

JET PRO^c



THE SPORICIDAL DISINFECTANT THAT KEEPS ON PROTECTING


JET PRO is a sporicidal chlorine dioxide foam. Each pull of the trigger targets surfaces where microbes are a risk. Not only that, an extra ingredient enables the foam to kill bacteria and enveloped viruses for up to 8 hours after application.

Disinfectants cannot be used every minute of every day. The combination of JET PRO's powerful chlorine dioxide chemistry with prolonged protection technology provides assurance for patients, staff and visitors, even after you have moved on to your next task.


PROLONGED PROTECTION FOR PATIENTS, STAFF & VISITORS



JET PRO is sporicidal in 1 minute.



Once dry, JET PRO places an invisible barrier on the surface, protecting the surface for up to 8 hours.



Touching the surface weakens the protective barrier, and regrowth of microorganisms will eventually occur.



JET PRO can then be re-applied to remove microorganisms and refresh the protective barrier.

WET KILL

DRY KILL

WET KILL

BOTH SWORD AND SHIELD TRISTEL ClO₂ + PROLONGED PROTECTION TECHNOLOGY

JET PRO is built upon the same chlorine dioxide (ClO₂) chemistry as the Tristel Medical Instrument Decontamination products. JET PRO incorporates two separate compartments that contain Base Solution (citric acid) and Activator Solution (sodium chlorite).

When mixed by pulling the JET PRO trigger, ClO₂ chemistry is generated. JET PRO contains an amphoteric surfactant for added cleaning performance. ClO₂ eliminates microorganisms by destroying their molecular structures. This means that harmful microbes cannot build resistance to ClO₂. An extra ingredient enables the foam to kill bacteria and enveloped viruses for up to 8 hours after application.

SPORICIDAL WIPES ON THE GO

Efficacy with the standard wipe used in EN 16615 means JET PRO can be paired with any dry wipe.

Use a microfibre cloth or a biodegradable dry wipe to create sporicidal wipes. You can stop using pre-wetted plastic wipes.

PROVEN EFFECTIVE IN 1 MINUTE, AGAINST:


SPORICIDAL	EN 17126
MYCOBACTERICIDAL	EN 14348
VIRUCIDAL	EN 14476
	PAS 2424*
YEASTICIDAL/FUNGICIDAL	EN 13624
	EN 13697
	EN 16615
BACTERICIDAL	EN 13727
	EN 16615
	PAS 2424

<i>Bacillus cereus</i>	<i>Aspergillus brasiliensis</i>
<i>Bacillus subtilis</i>	<i>Staphylococcus aureus</i>
<i>Clostridioides difficile</i>	<i>Enterococcus hirae</i>
Tuberculosis (TB) causing Mycobacteria	<i>Pseudomonas aeruginosa</i>
Poliovirus	<i>Escherichia coli</i>
Adenovirus	MRSA
Norovirus	CRE
SARS-CoV-2**	MDRAB
Bloodborne Viruses (e.g. HBV and HIV)***	VRE
<i>Candida auris</i>	
<i>Candida albicans</i>	

*Modified for virucidal efficacy against enveloped viruses.

READY TO USE APPLY, WIPE & WAIT

Each bottle of JET PRO delivers 570 doses of sporicidal disinfectant foam. **Simply apply, wipe and wait!**



ALSO AVAILABLE



JET

JET is a sporicidal chlorine dioxide foam. Each pull of the trigger targets surfaces where microbes are a risk.



FUSE

Designed for use on floors, FUSE combines cleaning action with the sporicidal power of chlorine dioxide.



TANK

TANK is a storage and distribution system fuelled by either Detergent or ClO₂ Sporicidal Disinfectant capsules.

ORDER NOW



**One Reusable Trigger
& Two JET PRO Bottles:
CCH011001**

Our team can offer on-site training, product rollout and audit support for your hospital. Get in touch to find out more.

cache[™]
Sporicidal Surface
Disinfection

+44 (0) 1638 721500 - mail@tristel.com
www.tristel.com

Copyright © Tristel Solutions
Mkt-Bro-1263-2
22/MAR/2022

Made in the UK by:

Tristel[™]
WE HAVE CHEMISTRY.

A representative sample of Tristel chlorine dioxide chemistry has been tested in accordance with EN 14476:2013+A2:2019, at a concentration of 20 parts per million (ppm). JET PRO has a chlorine dioxide concentration greater than 20ppm at point of use. *Based on virucidal efficacy validated in accordance with EN 14476:2013+A2:2019.